OptiBond™ **FL**

**Instruction For Use**

**Bonding Agent**

## **OptiBond™ FL**

### **Instructions for Use**

OptiBond FL is a two-part primer adhesive system indicated for direct light cured applications. OptiBond FL Adhesive is 48% filled with 0.6 microns barium glass, is fluoride releasing and is radiopaque. This unique technology provides the highest level of protection against microleakage, while sustaining high bond strengths to a variety of surfaces. The filled adhesive provides a true “structural bond” not found in unfilled adhesive systems. With OptiBond FL, adhesion promoters are carried in an ethanol solvent, diminishing the need for multiple coats.

OptiBond FL is available in both bottle and Unidose™ delivery systems.

Directions for using both delivery systems are as follows:

#### **Unidose Delivery System**

Unidose Delivery – For single patient use only, to prevent cross contamination between patients.

#### **To dispense OptiBond FL using the Unidose delivery system:**

- Open foil pouch along perforated edge – peel open the pouch.
- Remove the OptiBond FL Unidose yellow (FL Prime) and Unidose black (FL Adhesive) capsules and applicator brush from foil; discard foil.
- Grip Unidose capsule at wings.
- Twist each side in opposite directions until seal is broken completely.
- Discard side with the long stem (longer half).
- Dip applicator brush into opening, wetting tip completely.
- Discard Unidose device after use.

#### **Guidelines for bonding to enamel and/or dentin using the Unidose delivery system:**

1. Isolate prep and dry.
2. Place Kerr Gel Etchant\* with 37.5% phosphoric acid on dentin/enamel for 15 seconds. Rinse with water until etchant has been completely removed – approximately 15 seconds. Gently air dry for a few seconds being careful not to desiccate dentin.

3. Place applicator tip into the yellow Unidose capsule (OptiBond FL Prime). Apply material to the prepared enamel/dentin surfaces with a light scrubbing motion for 15 seconds. Gently air dry for approximately 5 seconds. At this point the dentin/enamel surface should have a slightly shiny appearance. Do not allow salivary contamination.  
**Note:** When restoring multiple teeth, use OptiBond FL Prime on all prepared surfaces before applying the Adhesive. If the applicator tip has been placed in the Adhesive DO NOT PLACE back into primer chamber - use a new applicator brush.
4. Using the same applicator brush, place tip into the black Unidose capsule (OptiBond FL Adhesive). Apply material to the prepared enamel/dentin surfaces with a light scrubbing motion for 15 seconds. Blow to margin or to thin if necessary using a light application of air.
5. Light cure \*\*.
6. Proceed with composite placement.

### **BOTTLE Delivery System**

#### **To dispense OptiBond FL Bottle delivery:**

- Dispense one drop of OptiBond FL Prime into disposable mixing well.
- Dispense one drop of OptiBond FL Adhesive into disposable mixing well.
- Use applicator brush to apply primer first and adhesive second.

#### **Guidelines for bonding to enamel and/or dentin using the Bottle delivery system:**

1. Isolate prep and dry.
2. Place Kerr Gel Etchant\* with 37.5% phosphoric acid on enamel and dentin for 15 seconds. Rinse with water until etchant has been completely removed – approximately 15 seconds. Gently air dry for a few seconds being careful not to desiccate dentin.
3. Dispense OptiBond FL Prime (Bottle #1) into disposable mixing well. Apply material over enamel and dentin surfaces with a light scrubbing motion for 15 seconds. Gently air dry for approximately 5 seconds. At this point the dentin surface should have a slightly shiny appearance. Do not allow salivary contamination.

4. Dispense and apply OptiBond FL Adhesive (Bottle #2) over enamel and dentin uniformly creating a thin coating. Blow to margin or to thin if necessary using a light application of air.
5. Light cure \*\*.
6. Proceed with composite placement.

Store at ambient temperature.

**\*CAUTION:** Contains phosphoric acid. Avoid contact with skin, eyes, and soft tissue. In case of contact with skin or eyes, flush immediately with water. Get medical attention for eyes. Do not take internally

**\*\*Recommended Cure Times when using Kerr Lights:**

Demi/Demi Plus, 5 seconds; L.E.Demetron II, 5 seconds; L.E.Demetron I, 10 seconds; or Optilux 501 in Boost mode, 10 seconds, Ramp Mode, 20 seconds, or Regular Mode, 20 seconds. For all other lights, see manufacturer's recommendation.

**CAUTION:** Uncured methacrylate resin may cause contact dermatitis and damage the pulp. Avoid contact with skin, eyes and soft tissue. Wash thoroughly with water after contact.

Kerr, OptiBond, and Unidose are trademarks owned by Kerr Corporation. Kerr, OptiBond, and Unidose are registered with the U.S. PTO and other countries.

**Limited Warranty - Limitation of Kerr's Liability**

Kerr's technical advice, whether verbal or in writing, is designed to assist dentists in using Kerr's product. The dentist assumes all risk and liability for damages arising out of the improper use of Kerr's product. In the event of a defect in material or workmanship, Kerr's liability is limited, at Kerr's option, to replacement of the defective product or part thereof, or reimbursement of the actual cost of the defective product. In order to take advantage of this limited warranty, the defective product must be returned to Kerr. In no event shall Kerr be liable for any indirect, incidental, or consequential damages.

EXCEPT AS EXPRESSLY PROVIDED ABOVE, THERE ARE NO WARRANTIES, BY KERR, EXPRESS OR IMPLIED, INCLUDING WARRANTIES WITH RESPECT TO DESCRIPTION, QUALITY, OR FITNESS FOR A PARTICULAR PURPOSE.

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**Kerr Australia Pty. Ltd.**
**Sybron Dental Specialties Japan, Inc.**