

Demetron Accessories



20879
Light Guide Holder
Holds six accessory light guides.



21042
Optics Maintenance Kit
After numerous autoclave cycles, brown discoloration may accumulate on polished ends of the light guide. Avoid by polishing ends every 50 cycles. Do not use to remove composite.



Your practice is *our* inspiration.™



Demetron
Light Guides

kerrdental.com 800.KERR.123

0609 5k

Demi only



921677
2mm

Curved
For layer curing of posterior estorations, inside gingival box preparation, and anterior tacking of laminate veneers.

Demi only



921676
4mm

Turbo Curved
For preventing overlapping on the gingiva in Class V.



20945
2mm

Curved
For layer curing of posterior restorations, inside gingival box preparation, and anterior tacking of laminate veneers.



21185
4mm

Turbo Curved
For preventing overlapping on the gingiva in Class V.



921535
8mm

Curved Amber
For general and posterior curing.



21020
8mm

Curved Turbo
Increases efficiency of light output and reduces overall curing time.



20941
8mm

Curved
For improved access in general and posterior curing.



952213
8mm

Turbo+
Increases efficiency of light output and reduces overall curing time.



902495
8mm

Turbo Extended
Extended length improves access to the posterior. Fits all Demetron curing lights.



20898
11mm

Curved
For improved access in general and posterior curing.



921551
11mm

Curved Turbo
Increases efficiency of light output and reduces overall curing time.



20812
13mm

Curved
For treating entire occlusal surface of molars.