

Kerr Jan/Feb Promotions 2012

January 1st - February 29th 2012



SonicFill™

new sonic-activated, single-step bulk fill composite



SonicFill is a unique bulk fill system comprised of a specially designed handpiece and new composite material in Unidose tips. SonicFill's composite incorporates a highly-filled proprietary resin with special modifiers that react to sonic energy.


Simply stated, there is no other technology like it on the market today! As sonic energy is applied through the handpiece, the modifier causes the viscosity to drop (up to 87%), increasing the flowability of the composite enabling quick placement and precise adaptation to the cavity walls. When the sonic energy is stopped, the composite returns to a more viscous, non-slumping state that is perfect for carving and contouring.

SonicFill Advantages


- Effortless placement
- Superior adaptation
- Easy contouring without stickiness
- No need for 2 mm incremental layering
- No final capping layer required
- 20 second full depth of cure

Application Time Savings

SonicFill vs Bulk Fill Flowables and Traditional Composite



Traditional layering technique
Several Layers



Bulk fill flowable
with universal cap
Two Layers



SonicFill System
ONLY ONE LAYER

SPECIAL OFFERS:

Buy SonicFill Intro Kit, **GET** an over the shoulder video by Dr. Ron Jackson and a Getting Started Kit **FREE!**

Buy 3 x Sonic Fill Composite Refills, **GET** 1 x Optibond XTR or 1 x Optibond Solo Plus, Bottle or Unidose **FREE!**

Your Practice is our inspiration™

Herculite Ultra Universal Nanohybrid Dental Composite




- **Optimal handling, non-slumping.** Superb ease of placement.
- **Long-lasting polishability.** Maintains lifelike luster long after initial restoration.
- **Herculite gold standard technology.** Proven 25+ years of experience.
- **High wear strength and mechanical properties.** Long-lasting and durable.
- **Chameleon quality.** Superior tooth-blending effect enabled by optimized filler size.
- **Now available in 21 shades.**

Premise Universal Nanofilled Composite



- **Ultra-Low Shrinkage.** With its 84% filler load, shrinkage is minimized to a low of 1.6% - up to 36% less than other brands.
- **Easy, Sustained Polish.** Superior luster retention over time.
- **Excellent Handling.** Highly sculptable and non-sticky.
- **Universal Application.** Strength for posteriors. Beauty for anteriors.

Premise Flowable



- **Clinically proven.** Reliable strength. Patented trimodal formula delivers excellent mechanical strength.
- **Ultra-low shrinkage.** 20% less shrinkage than leading brands, minimizing stress, sensitivity and microleakage.
- **Highly radiopaque.** Over 300% radiopacity - the highest of any flowable - makes X-ray detection easy.
- **Optimized viscosity.** Flows easily to fill cavities, yet holds its shape.
- **Gloss retention.** Polishability is easily achieved and sustained over time.

SPECIAL OFFERS:

Buy 3 x Herculite® Ultra or Premise™ or Premise Flowable GET 1 x Optibond XTR or 1 x Optibond Solo Plus (Bottle or Unidose) FREE!
Buy 5 x Herculite® Ultra OR Premise™ or Premise GET 2 of the same (Mix and Match OK. Unidose® or Syringe) FREE!*

Your Practice is our inspiration™

Vertise Flow



- **Self-adhering.** No separate bonding step.
- **Greatly reduces the chance of post-op sensitivity.** Shares the same characteristic inherent in self-etch materials.
- **High bond strengths to dentin and enamel.** Incorporates clinically proven OptiBond adhesive technology for excellent bonds.

Buy a Vertise Flow Intro Kit, Get a Trial Kit FREE!*

Vertise™ Flow, the first self-adhering composite powered by OptiBond®, greatly simplifies your direct restorative procedures by incorporating a bonding agent into a flowable. So, no need to bond separately. Fewer steps saves you time. And Vertise Flow shares the same characteristic inherent in self-etch materials: it reduces the chance of post-op sensitivity. Our leadership in resin restorative dentistry comes after decades of innovation in composites and expertise in adhesives - both of which come together in Vertise Flow. Visit www.kerrdental.com for the scientific research behind the technology.

Revolution Formula 2 Flowable Light Cure Composite



Revolution Formula 2™ is a flowable, light-cure, hybrid resin restorative. Building on the success of Revolution (the first flowable composite material on the market), Formula 2 has improved physical properties, fluoride release, more exact Vita®* shading and an expanded shade range.

SPECIAL OFFERS:

Buy 3 x Revolution® 2™ or Vertise Flow GET 1 x Optibond XTR or 1 x Optibond Solo Plus (Bottle or Unidose) FREE!
Buy 5 x Revolution® 2™ or Vertise Flow GET 2 of the same FREE!*

Your Practice is our inspiration™

Take 1 Advanced VPS Impression Material



- **Outstanding tear strength and elongation.** Strong, flexible material stretches around undercuts without tearing or distorting.
- **Hydrophilic wash material.** Displaces intraoral fluids to produce impressions with great detail and no voids or gaps.
- **Excellent dimensional stability.** Impression will not distort due to dimensional change or temporary deformation.
- **Widest selection of viscosities, set times and delivery methods.** Helps overcome common difficulties in impression taking, and matches doctor preference.

NEW - Rigid Viscosity - specifically for dual-arch impression trays!

Take 1® Advanced™ vinylpolysiloxane (VPS) impression material employs Kerr's long history of expertise in developing high-performance impression materials, along with recent technological advancements, to deliver a VPS material with the optimum combination of physical properties: strength, elasticity, dimensional stability and the ability to register detail in any environment. Furthermore, Take 1 Advanced comes in the widest selection of set times, delivery systems and consistencies available on the market today - so you always have the perfect combination of mechanical properties in a material that feels like it was designed just for your technique preferences.

SPECIAL OFFERS:

Buy 4 x Take 1® Advanced™ OR AlgiNot™ and AlgiNot FS™, GET 1 of the same Plus 1 x 15ml VPS Adhesive FREE!

*(Free! Good of Equal or Lesser Value. 50 ml Cartridges, Jars, Unidose® or 380 ml Volume Cartridges)

Take 1 Advanced Bite Registration



- **Highly stackable and no slump,** yet easy to extrude – stays where you put it, easy on the hand
- **Easy to trim,** no flaking, low deformation
- **Dimensionally stable**
- **Comparatively low contact angle** – more moisture tolerant to accurately record details
- **Work time 30 seconds,** total time 55 seconds – no time to move or distort, greater patient comfort
- **Durable, strong and accurate**
- **Fluffy and mousse-like** – easy to bite (patient comfort)

SPECIAL OFFER:

Buy 2 x Take 1® Advanced™ Bite Registration GET 1 x FREE!

Your Practice is our inspiration™

Maxcem Elite Self-Etch, Self-Adhesive Resin Cement



- **Bond strengths between 22-36 MPa.** Depending on substrates. High bond strengths without an adhesive.
- **Ease of use.** No need to refrigerate, easy cleanup, and no hand mixing required.
- **Compatibility with common substrates.** Anterior or posterior ceramics, PFMs and metal restorations, and CAD/CAM materials
- **Efficient self-cure mechanism.** Efficient self-cure mechanism. Sets up completely in the absence of light.

SPECIAL OFFERS:

Buy 2 x Maxcem Elite™ Refills **GET** 1 x Tempbond or Tempbond NE or Tempbond Clear (Tubes, Unidose or Auto Syringe) **FREE!**

NEW (MAXCEM ELITE BULK)!

Buy 2 x Maxcem Elite Bulk **GET** 1 x Maxcem Elite Bulk **FREE!**

NX3 Nexus Third Generation



- **Simplified delivery.** Dual-cure automix syringe eliminates hand mixing.
- **Light-cure applications.** Cement for veneers and indications requiring unlimited work time.
- **Bonds to all substrates.** Excellent adhesion to dentin, enamel, CAD/CAM blocks, ceramic, porcelain, resin and metal.
- **Self-etch or total-etch.** Bonding protocol compatibility. No dual-cure activator required.
- **Superior color stability.** Long-term esthetics for both dual-cure and light-cure cements.
- **Optimal handling.** Easy cleanup in gel state

SPECIAL OFFERS:

Buy 2 x NX3 Nexus™ Third Generation Refills - NX3 Dual-Cure Refills or NX3 Light-Cure Refills (Excludes Try-In Gels) (Mix and Match OK) **GET** 1 Tempbond or Tempbond NE or Tempbond Clear (Tubes, Unidose or Auto Syringe) **FREE!**

Your Practice is our inspiration™

OptiBond XTR Self-Etch, Light-Cure Universal Adhesive



- **Outstanding bond strengths to dentin and enamel.**
- **Self-etch.** Minimizes post-operative sensitivity, maximizes patient comfort.
- **Direct and indirect restorations.** A true universal adhesive – eliminates the need for multiple bonding agents.
- **Proven OptiBond™ Technology.** Dependable clinical performance and reliability.

OptiBond™ XTR self-etch, light-cure universal adhesive stands apart as the only bonding agent that embodies both the power and durability of a total-etch adhesive and the simplified protocol of a self-etch technique. Universal compatibility enables use for all direct and indirect restorations. The unique chemistry of an enhanced self-etch primer and optimized adhesive in this 2-bottle system brings outstanding bond strengths to uncut enamel and dentin.

SPECIAL OFFERS:

Buy 2 x Optibond XTR GET 1 x Herculite® Ultra OR Premise™ or Premise Flowable or Revolution® 2™ or Vertise Flow FREE!

Buy 5 x Optibond XTR GET 2 x Optibond XTR FREE!

OptiBond Solo Plus Single Component Total-Etch Adhesive



- **Clinically proven.** Universally accepted, industry-standard products. Proven performance worldwide.
- **Strong, durable bond.** 15%-filled material penetrates and reinforces dentin tubules to provide a long-lasting bond.
- **Delivery options.** Choose 5 ml bottle or convenient Unidose® delivery

OptiBond® Solo Plus™ is a single-component dental adhesive designed for both direct and indirect bonding applications. It is 15% filled with the same 0.4 micron filler found in the Kerr Point 4 composite. The filler not only reinforces the hybrid zone but also penetrates the dentin tubules, creating a true "structured bond" not found in unfilled or even "nano" filled adhesive systems. This unique technology provides the highest level of protection against microleakage, while sustaining high bond strengths to a variety of surfaces.

SPECIAL OFFERS:

Buy 2 x Optibond Solo Plus or Optibond XTR or Optibond AIO

GET 1 x Herculite® Ultra or Premise™ or Premise Flowable or Revolution® 2™ or Vertise Flow FREE!

Buy 5 x Optibond Solo Plus or Optibond XTR or Optibond AIO

GET 2 x Optibond Solo Plus or Optibond XTR or Optibond AIO FREE!

Your Practice is our inspiration™

Blu/White Diamonds (NOW AVAILABLE IN 5 PACKS)



Unique diamond application process ensures excellent diamond particle retention and cutting surface throughout multiple uses

- A wide selection of grit sizes, including Extra Coarse for rapid gross reduction
- Wide selection of shapes to accommodate any clinical indication or prep design
- Outstanding cutting efficiency reduces chatter and instrument wear and tear

BluWhite® Diamond Burs cut faster and clean easier than other leading manufacturers due to the unique grit application process, which helps reduce clogging and improves cutting efficiency. With a wide selection of shapes and grits, including unique extra coarse diamond particles (40% larger than competitors), BluWhite diamonds offer a great option for any operative or prosthodontic procedure time after time.

SPECIAL OFFER:

Buy 10 x (5 Packs) of Blue White Diamond Burs **GET** 10 x (5 Packs) **FREE!**

Kerr/Jet Beaver Carbide Burs



New Single-piece construction of our Solid Carbide burs improves strength and reliability. Simultaneously increases corrosion resistance.



The Quality Name in Carbides

- **Dental Excavating Burs:** designed for removal of carious tooth structure and shaping of the cavity in preparation for restorative materials
- **Trimming & Finishing Burs:** designed for the shaping of restorative materials to tooth Anatomy and desired Aesthetic pearance.
- **Fine Finishing Burs:** designed for the final Finish on Restorations to achieve a surface finish acceptable Visually and Aesthetically.
- **Fine Cross Cut Burs:** designed originally for the removal of soft and hard filling materials and crown removal. Faster cutting, reduced chatter, more control
- **Oral Surgical Burs:** designed with an extra shank length for oral surgery and difficult, hard to reach areas.

Kerr/Jet Beaver Carbide Burs Available in five packs or 100 packs

SPECIAL OFFERS:

Buy 10 x (5 Packs) Kerr Jet Carbide Burs, Dental Excavating Burs (6 Blades) Trimming & Finishing (12 Blades) & Fine Cross Cut (30 blades)
GET 10 x (5 Packs) **FREE!**

Buy 4 x (100 Packs) of Jet Carbide Burs **GET** 1 x (100 Pack) **FREE!**

Your Practice is our inspiration™

Kerr Demi Plus LED Curing Light



- **Advanced Battery Design.** The positive contact charger guarantees a snug fit between the contacts on the handpiece and the charger base allowing you to get the most out of every charge.
- **Enhanced Body Design.** The exterior is made from Valox resin, a high-performance plastic that will protect and extend the life of your unit.
- **Extended Light Guide.** The Extended Turbo Light Guide enables better application with improved access (+28% additional reach) for posterior procedures.
- **Longer Warranty.** Our 3year warranty ensures that you get the most value out of your premium investment.

SPECIAL OFFERS:

Buy 2 x Demi™ Plus LED Curing Light, **GET 1 x FREE!**

Buy 1 x Demi™ Plus LED Curing Light, **GET 1 x Handpiece and Battery FREE!**

Ultra Waxter 2



The Ultimate electronic waxer just got better, with dual spatulas for increased productivity, convenient tip holder for eight tips. Easy to use menu function, such as quick heat and programmable temperature presets.


Ordering Information: 34844 Ultra-Waxer 2
 Contains: Ultra-Tips #35167 & #35168,
 2 x Ultra-Spatula Cord/handle # 35177,
 Power Cord, Instructions



Tip Order #	Tip Description
KE-35167	Small PKT, Blue
KE-35168	Large PKT, Teal
KE-35169	Beavertail, Purple
KE-35170	Small # 7, Grey
KE-35171	Large # 7, Green
KE-35172	Needle, Red
KE-35173	1/2 Hollenback, Black
KE-35174	Small Denture, Pink
KE-35175	Lge Denture Spoon, Or
KE-35176	Denture Spoon, Yellow
KE-35177	Ultra-Spatula Cord/Handle

SPECIAL OFFER: Buy 1 x Ultra-Waxer 2, **GET 2 x Ultra Tips FREE!**

Galetti Articulator



The world's finest plasterless, fixed plate articulator. Made in Italy, the Galetti Articulator is all metal constructed. It provides for easy, speedy firm grasp of models of any size, controlling the correct occlusion with maximum precision. Neither plaster work nor accessories are required. Easy to clean.

Order Information
 32685

SPECIAL PROMOTIONAL PRICE \$880⁰⁰

Your Practice is our inspiration™

Opti 1 Step Polisher

Enhanced polishing results for composite restorations: Smooth surface and high gloss, in just one step



- Increases surface gloss
- Eliminates scratches on composite restorations
- Reduces surface roughness and bacterial retention
- Delivers comparable results to multi-step systems and procedures
- Perfect control of abrasiveness on all surfaces due to optimised flexibility
- Allows a variable abrasion by changing pressure
- Gold plated mandrel and unique abrasive matrix for greater durability
- Autoclavable, reusable
- More cost efficient compared to traditional procedures

Order codes:

- 8000 Opti 1 Step Polisher Assortment Kit
- 8001 Opti 1 Step Polisher Refill Flame
- 8002 Opti 1 Step Polisher Refill Cup
- 8003 Opti 1 Step Polisher Refill disc

SPECIAL OFFERS:

Buy 1 x Opti1 Step Polisher Assorted Kit, **GET** 1 x 3 Pack of Optishine Polishers **FREE!**

Buy 2 x Opti1 Step refill Kits, **GET** 1 x Opti1 Step Refill **FREE!**

OptiDisc®



The high performing disc, now even easier to choose!

- OptiDisc® finishing and polishing system gives your restorations a final polish equal to that of natural dentition.
- Kerr OptiDisc® system is the perfect combination of different grit size discs for contouring, finishing and high gloss polishing for composite, glass ionomers, amalgams, semiprecious and precious metals.

SPECIAL OFFERS:

Buy 3 x Optidisc Refills, **GET 1 x FREE!**

Buy 1 x Optidisc Kit (Item # 4200), **GET 1 x Optidisc Refill or 1 x 3 Pack of Optishine Polishers FREE!**

(Sample Packs Available on Request (Limited to one per doctor))

MetaFix™ All-in-One Matrix System



The new All-in-One MetaFix™ matrix, with the innovative integrated tightening/opening system, is the perfect solution for Class II – MO/OD/MOD composite fillings in the posterior area.

Order codes:

- Item 3600 Metafix Assorted Kit
- Item 3601 MetaFix Refill – Small (50 pcs)
- Item 3602 MetaFix Refill – Medium (50 pcs)
- Item 3603 MetaFix Refill – Large (50 pcs)




SPECIAL OFFER:

Buy 1 x Metafix Assort Kit, **GET 1 x Metafix Refill FREE!**

Your Practice is our inspiration™

X-Ray Film Holders



X-Ray Film Holder System/X-Ray Film Holders

- **Exact positioning** of the film or phosphor plate relative to the tooth – no bending and thus no distortion
- **No ocular guesswork** – exclusive, automatic reliance on the extraoral aiming and centring device
- **Reproducible images** (quasi standardisation)
- **Optimally clear radiographs**
- **Constant X-ray quality over time** – whoever takes the radiographs
- **Simplicity of the system in use** – saving time


SPECIAL OFFER:

Buy 1 x X-Ray Film and Phosphor Plate Holders Testset Item #1720

Contents: 6 X-ray holders with ring (1 Kwik-Bite, 1 Super-Bite Anterior, 1 Super-Bite Posterior, 1 Paro-Bite, 1 Endo-Bite Anterior, 1 Endo-Bite Posterior), 6 Centring Devices

GET 1 x Packet of Optdisc Mini Kit RR Value \$100.00 **FREE!**

Kwik Bite® Senso & Super Bite® Senso



Easy in. Easy out. Perfect images no matter which sensor you use

Kwik-Bite® & Super-Bite® Senso are universal sensor holders for direct digital imaging of horizontal bitewings and periapical radiographs. Thanks to the unique clamping device Senso holders are compatible with all commonly-used brands, types and sizes of sensors. Autoclavable at 134°C, at least 3 min.

SPECIAL OFFER:

Buy 1 x Testset Sensor Holder System Item # 2720

Contents: 5 sensor holders with ring (1 Kwik-Bite Senso, 1 Super-Bite Senso Anterior, 1 Super-Bite Senso Posterior, 1 Endo-Bite Senso Anterior, 1 Endo-Bite Senso Posterior), 5 Centring Devices

GET 1 x Box of SENSOR COVER 3.5 X 20.5cm(500) Item #3735-1 or 1 x Box SENSOR COVER 4 X 21cm (500) **FREE!**

Your Practice is our inspiration™

Presenting

Dr Marc Geissberger

D.D.S., M.A., B.S, CPT

Professor and Chairman of the Department of Restorative Dentistry at the University of the Pacific, School of Dentistry in San Francisco, California.

7
CPD
HOURS



Cutting Edge Techniques and Materials to Maximize Success of Anterior & Posterior Composite Resin Restorations

Hands-On Course Outline:

The use of composite resin material as a direct restorative for anterior and posterior restorations has dramatically increased in the market place. With self-etching technology emerging as a viable alternative to traditional techniques and the development of new technology to assist clinicians with bulk fill techniques, it is easy to get confused and fall behind. What materials will you advocate for use in your practice? What has the literature been revealing on materials in current use? Should practitioners change their current techniques? What preparation design will you utilize?

This lecture/hands-on course will attempt to make sense of all of the products and techniques available in today's market place and will focus on materials and their applications as well as a brief review of some of the literature. Emphasis will be placed on new developments in the composite area as well as outline what developments to watch for in the coming years.

During this program practitioners will learn:

- how to utilize self etching principles and how they are applied to composite resins
- how to reproduce anatomical form with every restoration
- how to successfully and predictably place a bulk-filled restoration
- proper placement techniques that maximize clinical outcomes and esthetics
- how to select and place an appropriate matrix system how to create and maintain correct posterior anatomy and contours
- how to finish composite restorations to create life-like results

Full Day Lecture Outline

The use of composite resin material as a direct restorative for anterior and posterior restorations has dramatically increased in the market place. With self-etching technology emerging as a viable alternative to traditional techniques and the development of new technology to assist clinicians with bulk fill techniques, it is easy to get confused and fall behind.

What materials will you advocate for use in your practice? What has the literature been revealing on materials in current use? Should practitioners change their current techniques? What preparation design will you utilize?

This lecture course will attempt to make sense of all of the products and techniques available in today's market place and will focus on materials and their applications as well as a brief review of some of the literature. Emphasis will be placed on new developments in the composite area as well as outline what developments to watch for in the coming years.

During this program practitioners will understand:

- self etching principles and how they are applied to composite resins
- what is required of materials to be utilized for bulk-filled and the differences between several current materials
- how to successfully and predictably place a bulk-filled restoration
- proper placement techniques that maximize clinical outcomes and esthetics
- how to select and place an appropriate matrix system
- how to create and maintain correct posterior anatomy and contours

Registration Form Yes, I would like to attend:

Hands-On Course Dates:

- Melbourne** Monday Feb 27th 2012 Number attending
Brisbane Wednesday March 7th 2012 Number attending

Lecture Dates:

- Auckland** Friday Feb 24th 2012 Number attending
Melbourne Wednesday Feb 29th 2012 Number attending
Adelaide Monday March 5th 2012 Number attending
Cairns Friday March 9th 2012 Number attending

Registration Time: 8.00 am to 8.30 am

Full Day Hands-On Course (Restricted to 30):

Includes Half/day lecture & Half/day hands-on Workshop. Material and instruments

Cost: \$660.00 inclusive of GST **Early Bird Price \$550.00** Inc of GST

Register Before February 1st 2012 to receive the Early Bird price
 (Kerr Composite & Bond Kit Included in price Valued at \$250.00)

Full Day Lecture Program with demo

Cost: \$396.00 inclusive of GST **Early Bird Price \$297.00** Inc of GST

Register Before February 1st 2012 to receive the Early Bird price

Please post or fax completed form to:

KERR AUSTRALIA Unit 10, 112-118 Talavera Rd, North Ryde, NSW 2113
Tollfree: 1800 643 603 Sydney Ph: (02) 8870 3000 • Fax: (02) 9870 7600

Name: _____

Address: _____

Postcode: _____ Phone: _____

Email: _____

Payment Details: Amount Enclosed: _____

Cheque Enclosed Visa MasterCard Diners Club American Express

Expiry Date: [] [] [] [] [] [] [] []

Holder's Name: _____

Signature: _____

Kerr Jan/Feb 2012 Redemption Coupon

Impression Materials

- Buy** 4 x Take 1 Advanced or AlgiNot™ or AlgiNot FS™
Get 1 of the same Plus 1 x 15ml VPS Adhesive FREE
 Cartridges, Jars, Unidose, 380 ml Volume
 Wash, Medium, Tray, Rigid Tray
 Super Fast/Set Fast /Set Reg/Set
- Buy** 2 x Take 1 Advanced Bite Reg, **Get** 1 x FREE

Composite

- Buy** 3 x Sonic Fill Composite Refills
Get 1 x Optibond XTR or 1 x Optibond Solo Plus FREE
 Bottle / Unidose
- Buy** 1 x SonicFill Intro Kit
Get an over the shoulder video by Dr. Ron Jackson and a Getting Started Kit FREE
- Buy** 3 x Herculite Ultra, Premise, Premise Flowable, Revolution 2 or Vertise Flowable
Get 1 x Optibond XTR or 1 x Optibond Solo Plus FREE
 Unidose / Syringe
- Buy** 5 x Herculite Ultra, Premise, Premise Flowable, Revolution 2 or Vertise Flowable
Get 2 of the same FREE
- | Qty | Brand | Shade | <input type="checkbox"/> Unidose / <input type="checkbox"/> Syringe |
|-------|-------|-------|---|
| _____ | _____ | _____ | <input type="checkbox"/> Unidose / <input type="checkbox"/> Syringe |
| _____ | _____ | _____ | <input type="checkbox"/> Unidose / <input type="checkbox"/> Syringe |
| _____ | _____ | _____ | <input type="checkbox"/> Unidose / <input type="checkbox"/> Syringe |

Cements

- Buy** 2 x Maxcem Elite™ Refills, or NX3 Dual-Cure Refills or NX3 Light-Cure Refills (Excludes Try-In Gels)
Get 1 x **Tempbond** Unidose, Automix, Tubes or **Tempbond NE** Unidose, Automix, Tubes or **Tempbond Clear**
- Buy** 2 x Maxcem Elite Bulk **Get** 1 x Maxcem Elite Bulk FREE

Bonds

- Buy** 2 x Optibond Solo Plus, Optibond XTR, Optibond AIO,
Get 1 x Herculite® Ultra or Premise or Premise Flowable or Vertise Flow FREE Bottle / Unidose
- Buy** 5 x Optibond Solo Plus, Optibond FL, Optibond AIO,
Get 2 x Optibond Solo Plus, Optibond FL, Optibond AIO, FREE
 Unidose Bottle

Small Equipment

- Buy** 2 x Demi Plus LED Curing Light **GET** 1 x FREE
- Buy** 1 x Demi Plus LED Curing Light
GET 1 x Handpiece and Battery FREE
- Buy** 1 x Ultra-Waxer 2 **GET** 2 x Ultra Tips FREE
 Please indicate: Tip # _____ Tip # _____

Kerr/Hawe Products

- Buy** 3 x Optidisc Refills **Get** 1 x FREE (Mix and Match OK) Item # _____
- Buy** 1 x Optidisc Kit
Get 1 x Optidisc Refill Item #4200 or 1 x 3 Pack of Optishine Polishers FREE
- Buy** 1 x Opti 1 Step Polisher Assorted Kit
Get 1 x Opti 1 Step Refill Cup Disc Point FREE
- Buy** 1 x Metafix Assort Kit **Get** 1 x Metafix Refill Item # _____ FREE
- Buy** 1 x Testset Sensor Holder System Item # 2720
Get 1 x SENSOR COVER 3.5 X 20.5CM(500)Item #3735-1 or 1 x Box SENSOR COVER 4 X 21CM (500) Item #3736-1 FREE
- Buy** 1 x X-Ray Film Holders Testset Item #1720
Get 1 x Packet of Optdisc Mini Kit FREE

Kerr Burs

- Buy** 10 x (5 Packs) Blu/White Burs **Get** 10 x (5 Packs) FREE #s _____
- Buy** 10 x (5 Packs) Jet Cardide Burs **Get** 10 x (5 Packs) FREE #s _____
- Buy** 4 x (100 Packs) of Jet Carbide Burs
Get 1 x (100 Pack) FREE #s _____

Kerr August 2011 Promotions Only Available in Australia/New Zealand

Please print **PROMO CODE N2027-1** legibly on your invoice.
 To process your free good(s), fax your coupon with proof of purchase to (02) 9870 7600 or email and attach a PDF of your proof of purchase to kerraustpromos@sybrondental.com

Name: _____

Office Name: _____

Address: _____

Postcode: _____

Phone: _____

Email: _____

I would like to receive promotional, technical, and product information via email. Please send me an email so that I can opt-in to Kerr eNews.*

*This information is not sent out in an unsolicited manner and is only sent to users that request the information. All email newsletters provide information on how to unsubscribe from email distribution lists at any time. The required purchase must be reflected on one preferred Kerr dealer invoice dated between January 1st - February 29th, 2012. Coupons and invoices received after March 15th, 2012 will not be honoured. Multiple redemptions on a single invoice are valid, but only one redemption of each offer outlined herein is allowed per office for non-dental goods. Offer cannot be combined with any other Kerr or dealer special. Allow 4 to 6 weeks for delivery of free goods. Free goods available while supplies last. Incomplete or illegible coupons or invoices will not be honoured. Kerr reserves the right to discontinue or change this offer at any time.

Your Practice is our inspiration™